




HOOD EXHAUST RATINGS:
 7.92 SQ FT OPENING AT FULL OPEN @ 100 LFM, 792 CFM
 4.75 SQ FT OPENING AT 18" OPEN @ 100 LFM, 475 CFM
 .33 STATIC PRESSURE AT FUME HOOD DUCT COLLAR



QTY ONE- AIR FOIL FUME HOOD REQ'D AS SHOWN



LONGO
 Inc.

100 HILLTOP ROAD
 RAMSEY, NJ 07446
 PHONE: 201.825.1500
 FAX: 201.825.4784
 WWW.LONGOLABS.COM

SHEET DESCRIPTION:

DETAIL

DATE: 2005.8.4

PROJECT #:

00000

SCALE: 1/2" = 1'-0"

DESIGNED BY: SEAN ONDICH

CHECKED BY: ANTHONY STELLATOS

ANTHONY@LONGOINC.COM

SHEET #:

D-001