



**NOTES:**

- 1) EXTERIOR- COLD ROLL STEEL W/ACID RESISTANT POWDER COATED FINISH. COLOR TO BE \_\_\_\_\_
- 2) INTERIOR- WHITE POLYRESIN MATERIAL.
- 3) FIXTURES- DECK MOUNT- POLISHED CHROME- (1) CW GOOSENECK (WATERSAVER P/N #L614) & (1) GAS (WATERSAVER #L4200-132A) SUPPLIED BY AMS. (SHIP LOOSE FOR DECK MOUNT) (PRE-PIPING BY OTHERS)
- 4) ALL FINAL CONNECTIONS BY RESPECTIVE TRADES.
- 5) ELECTRICAL- (1) LT SW & (1) 115V DUP SUPPLIED BY AMS. (SHIPPED LOOSE FOR MOUNTING IN BASE BY OTHERS)
- 6) DISHED EPOXY RESIN COUNTERTOP W/ SINK SUPPLIED BY COLLEGE DALE
- 7) HOOD SPECS:  
 1250 CFM @ 60 LFM WITH ONE SASH OPEN  
 400 CFM @ 60 LFM WITH ONE SASH OPEN @ 18"  
 .45 STATIC PRESSURE AT HOOD DUCT COLLAR



QTY ONE- DEMONSTRATION HOOD FUME HOODS REQ'D AS SHOWN



DEALER: LONGO ASSOCIATES - 100 HILLTOP ROAD - RAMSEY, NEW JERSEY - 07446  
 PHONE: 201-825-1500 FAX: 201-825-4784 WWW.LONGOSCHOOLS.COM

DRN.	DATE		
CHK'D	DATE	ITEM	
SCALE		SHEET	OF
REVISIONS		DRAWING NO.	

CUSTOMER