




HORIZONTAL SECTION

HOOD EXHAUST RATINGS:
 7.92 SQ FT OPENING AT FULL OPEN @ 100 LFM, 792 CFM
 4.75 SQ FT OPENING AT 18" OPEN @ 100 LFM, 475 CFM
 .33 STATIC PRESSURE AT FUME HOOD DUCT COLLAR
 (BASED W/ ONE SASH OPEN)



SECTION THRU

QTY ONE--48-MDS-DOUBLE-SIDED FUME HOOD REQ'D AS SHOWN



LONGO
INC.

100 HILLTOP ROAD
 RAMSEY, NJ 07746
 PHONE: 201.825.1500
 FAX: 201.825.4784
 WWW.LONGOLABS.COM

SHEET RESPONSIBILITY
 DETAIL

DATE: 2005.8.4

PROJECT # 00000 SCALE 1/2" = 1'-0"

DESIGNED BY SEAN ONDICH

CHECKED BY ANTHONY STELLATOS

ANTHONY STELLATOS

SHEET # D-011