

6. 76230 TRIFACIAL SERVICE CENTER

(3) sides at 34" ea. and (3) sides at 20" ea.-55 1/2" at widest dimension. Countertop and sink are made of solid epoxy resin. Hexagonal drop in sink with raised rib basin is 30"x26"x6.9".

Center pedestal is 14 3/4" square, made of 14 gauge formed steel and provides for attachment of support structure. A 12"x25" access panel is located on one side for access to mechanical connections.

Support understructure is totally enclosed and made of molded fiberglass. Understructure is designed to mount various types of re-locatable storage units. Removable access panels are located directly below service connections. Molded fiberglass, when subjected to the ASTM D-635 test procedure, shall exhibit the following flammability characteristics:

ASTM D-635	NON-BURNING
HLT-15 (AVERAGE)	60
OXYGEN INDEX (AVERAGE)	30.0
SMOKE DEVELOPMENT (AVERAGE)	300

Notarized letter signifying manufacturer's products comply with flammability requirements furnished when required.

Service fixtures include (3) G.F.I. duplex electrical outlets and (3) solid cast brass vandal resistant combination hot water/cold water/gas fixtures. Electrical outlets must be mounted above adjoining table countertop height. Internal piping and wiring are not included.

When installed properly, unit should support around the perimeter a 250lb. load.