



Oak (standard)



Maple



Zebra Stripe Bamboo



Plastic Laminate

## NOT LIKE THE OTHERS

Our durable, laboratory-grade casework is perfect for science applications from kindergarten to college. Our exclusive thru-bolt construction has solid wood frame members that lock corners together, creating superior cabinet construction not matched by our competitors. Sturdy, adaptable and constructed from the highest grades of wood in the industry, Sheldon's casework can be counted on for years of durable and reliable wear. You can choose from our many contemporary door and drawer fronts to suit your styling needs.

### ***Sheldon Casework Features:***


- ✓ High Grade Interior **And** Exterior Wood
- ✓ Durable Thru-bolt Construction
- ✓ Full-framed Tops **And** Bottoms
- ✓ Vertical Match Grain Fronts
- ✓ Various Styles **And** Finishes Available

***Innovative • Reliable • On Time • Complete***

 **Sheldon**  
LABORATORY SYSTEMS  
since 1898

***Casework***

For detailed information about  
Sheldon products,  
call 1.800.531.7604 or visit  
[www.sheldonlabs.com](http://www.sheldonlabs.com)



**Wood Veneers** - Oak, maple and other wood species available.

**Hardware** - Powder-coated epoxy (black, chrome or white) or stainless steel. Additional pull styles available upon request.

**Finishes** - Variety of finishes to choose from; custom color matches available upon request.

**Drawer Slides** - Self closing, soft closing or full extension available upon request.

**Tops** - Epoxy resin available in a variety of standard colors. Black phenolic resin, plastic laminate tops in standard colors.

**Drawer Boxes** - Standard with chuck and bore construction and woodgrain edgbanding. Also available with dovetail construction.