



HOOD EXHAUST RATINGS:
 922 SQ FT OPENING AT FULL OPEN @ 100 LFM, 922 CFM
 553 SQ FT OPENING AT 18" OPEN @ 100 LFM, 553 CFM
 .30 STATIC PRESSURE AT FUME HOOD DUCT COLLAR



QTY 1 - THIN WALL FUME HOODS REQ'D AS SHOWN



LONGO
inc.

100 HILLTOP ROAD
 RAMSEY, NJ 07446
 PHONE: 201.825.1500
 FAX: 201.825.4784
 WWW.LONGOLABS.COM

SHEET DESCRIPTION	
DETAIL	
DATE	2005.8.4
TITLE	
PROJECT #	00000
SCALE	1/2" = 1'-0"
DRAWN BY	SEAN ONDICH
CHECKED BY	ANTHONY STELLATOS
SHEET #	D-005