



Assemblies

Island Labs



2 PERSON WORKSTATION - Single faced tables that provide a fully equipped science center for two people. These single faced units contain a solid epoxy resin sink along with two faucets. Every unit contains enough gas cocks for each person along with GFI AC duplex electrical receptacles, rod sockets and book compartments. The cabinet contains a cupboard with an adjustable shelf. 66"L x 30"W x 36"H.


ISL-210A-PL 2 Person Workstation



Door Option A

2 PERSON WORKSTATION WITH FLAT TOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

ISL-210ACT-PL 2 Person Workstation with Flat Top



4 PERSON WORKSTATION - Single faced tables that provide a fully equipped science center for four people. These single faced units contain two solid epoxy resin sinks along with four faucets. Every unit contains enough gas cocks for each person along with GFI AC duplex electrical receptacles, rod sockets and book compartments. The two cabinets each contain a cupboard with an adjustable shelf. 11'L x 30"W x 36"H.

ISL-410A-PL 4 Person Workstation


4 PERSON WORKSTATION WITH FLAT TOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

ISL-210ACT-PL 4 Person Workstation with Flat Top



2 PERSON WORKSTATION - Single faced tables that provide a fully equipped science center for two people. These single faced units contain a solid epoxy resin sink along with two faucets. Every unit contains enough gas cocks for each person along with GFI AC duplex electrical receptacles, rod sockets and book compartments. The cabinet contains 4 drawers and a cupboard. 66"L x 30"W x 36"H.


ISL-210B-PL 2 Person Workstation



Door Option B

2 PERSON WORKSTATION WITH FLAT TOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

ISL-210BCT-PL 2 Person Workstation with Flat Top



4 PERSON WORKSTATION - Single faced tables that provide a fully equipped science center for four people. These single faced units contain two solid epoxy resin sinks along with four faucets. Every unit contains enough gas cocks for each person along with GFI AC duplex electrical receptacles, rod sockets and book compartments. The two cabinets each contain 4 drawers and a cupboard. 11'L x 30"W x 36"H.


ISL-410B-PL 4 Person Workstation

4 PERSON WORKSTATION WITH FLAT TOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

ISL-410BCT-PL 4 Person Workstation with Flat Top


Assemblies

Island Labs



4 PERSON WORKSTATION - Double faced tables that provide a fully equipped science center for four people. These double faced units contain a solid epoxy resin sink along with two faucets. Every unit contains enough gas cocks for each person along with GFI AC duplex electrical receptacles, rod sockets and book compartments. The front and back of the cabinet each contain a cupboard with an adjustable shelf. 66"L x 48"W x 36"H.


ISL-420A-PL 4 Person Workstation



Door Option A

4 PERSON WORKSTATION WITH FLAT TOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

ISL-420ACT-PL 4 Person Workstation with Flat Top




8 PERSON WORKSTATION - Double faced tables that provide a fully equipped science center for eight people. These double faced units contain two solid epoxy resin sinks along with four faucets. Every unit contains enough gas cocks for each person along with GFI AC duplex electrical receptacles, rod sockets and book compartments. The front and back of the two cabinets each contain a cupboard with an adjustable shelf. 11'L x 48"W x 36"H.

ISL-820A-PL 8 Person Workstation


8 PERSON WORKSTATION WITH FLAT TOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

ISL-820ACT-PL 8 Person Workstation with Flat Top



4 PERSON WORKSTATION - Double faced tables that provide a fully equipped science center for four people. These double faced units contain a solid epoxy resin sink along with two faucets. Every unit contains enough gas cocks for each person along with GFI AC duplex electrical receptacles, rod sockets and book compartments. The front and back of the cabinet each contain 4 drawers and a cupboard. 66"L x 48"W x 36"H.


ISL-420B-PL 4 People Workstation



Door Option B

4 PERSON WORKSTATION WITH FLAT TOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

ISL-420BCT-PL 4 Person Workstation with Flat Top



8 PERSON WORKSTATION - Double faced tables that provide a fully equipped science center for eight people. These double faced units contain two solid epoxy resin sinks along with four faucets. Every unit contains enough gas cocks for each person along with GFI AC duplex electrical receptacles, rod sockets and book compartments. The front and back of the two cabinets each contain 4 drawers and a cupboard. 11'L x 48"W x 36"H.


ISL-820B-PL 8 Person Workstation

8 PERSON WORKSTATION WITH FLAT TOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

ISL-820BCT-PL 8 Person Workstation with Flat Top

Assemblies

Island Labs




12 PERSON WORKSTATION - Double faced tables that provide a fully equipped science center for 12 people. These double faced units contain three solid epoxy resin sinks along with six faucets. Every unit contains enough gas cocks for each person along with GFI AC duplex electrical receptacles, rod sockets and book compartments. The front

and back of the three cabinets each contain a cupboard with an adjustable shelf. 16½'L x 48"W x 36"H.

ISL-1210A-PL 12 Person Workstation

12 PEOPLE WORKSTATION WITH FLAT TOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

ISL-1210ACT-PL 12 People Workstation with Flat Top




12 PERSON WORKSTATION - Double faced tables that provide a fully equipped science center for 12 people. These double faced units contain three solid epoxy resin sinks along with six faucets. Every unit contains enough gas cocks for each person along with GFI AC duplex electrical receptacles, rod sockets and book compartments. The front

and back of the three cabinets each contain 4 drawers and a cupboard. 16½'L x 48"W x 36"H.

ISL-1220B-PL 12 Person Workstation

12 PERSON WORKSTATION WITH FLAT TOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

ISL-1220BCT-PL 12 Person Workstation with Flat Top




Door Option A



Door Option B

Assemblies


Lab Benches



DOUBLE SIDED ISLAND LAB CENTER WITH RAISED SERVICE TURRET - Excellent choice for the physics lab. Designed for eight people. Raised turret houses four single gas cocks, four duplex 110-V AC electrical outlets. 1" flat phenolic resin countertop. Raised turret is selected hardwood with 1" phenolic resin cap. Substitute any base cabinets totaling 12'- 0" long. Turret 14"W x 12"D x 12"H.

LRB-203-PL

W	12' - 2"
H	36 to countertop / 48 with raised turret
D	48



DOUBLE SIDED ISLAND LAB CENTER WITHOUT KNEE SPACE - Designed for eight to twelve people. 1" epoxy resin countertop with trough sink running entire length of the assembly, draining into end sink. Reagent shelving is 3/4" phenolic resin. Includes six cold water cocks, six single gas cocks, 6 duplex 110-V AC electrical outlets and one cold water bibb cock at end sink. Service piping and wiring not included. Substitute any base cabinets totaling 12' - 0" long. End sink 22"L x 13½"W x 12½"D.


LRB-218-PL

W	12' - 2" not including end sink / 11' - 11" with reagent shelf
H	36 to countertop / 54 with reagent shelf
D	54 countertop / 12 with reagent shelf

Combination options, sizes and lengths can be customized to fit your requirements.

Assemblies


Lab Benches



DOUBLE SIDED ISLAND LAB CENTER WITH KNEE SPACE - Designed for eight people. 1" flat epoxy resin countertop (no trough sink). Includes three A-way gas cocks, three double sided electrical pedestals with one duplex 110-V AC outlet each side, one cold water deck mounted gooseneck fixture and one end sink. Reagent shelf 3/4" phenolic resin. Substitute base cabinets with any 24" x 22"D base cabinet. End sink 22"L x 13½"W x 12½"D.

LRB-230-PL

W	16' - 0" not including end sink / 14' - 9" with reagent shelf
H	36 to countertop / 54 with reagent shelf
D	54 countertop / 12 with reagent shelf



DOUBLE SIDED ISLAND LAB CENTER WITH KNEE SPACE - Designed for eight people laboratory center. Each person has individual storage at his lab station. 1" epoxy resin countertop with trough sink running the entire length of the assembly draining into the end sink. Reagent shelving is 3/4" phenolic resin. Includes eight cold water cocks, eight single gas cocks, 8 duplex 110-V AC electrical outlets and one cold water bibb cock at end sink. Service piping and wiring in reagent rack not included. Substitute base cabinets with any 24"W x 22"D base cabinets. End sink 22"L x 13½"W x 12½"D.


LRB-240A-PL

W	16' - 0" not including end sink / 15' - 9" with reagent shelf
H	36 to countertop / 54 with reagent shelf
D	54 countertop / 12 with reagent shelf

Combination options, sizes and lengths can be customized to fit your requirements.

Assemblies

Lab Benches



DOUBLE SIDED ISLAND LAB CENTER WITHOUT KNEE SPACE - Designed for eight to twelve people. 1" flat epoxy resin countertop (no trough sink). Includes three 4-way gas cocks, three double sided electrical pedestals with one duplex 110-V AC outlet each side, one cold water deck mounted gooseneck fixture and one end sink. Reagent shelf 3/4" phenolic resin. Substitute any base cabinets totaling 12' - 0" long. End sink 22"L x 13½"W x 12½"D.


LRB-450-PL

W	12' - 2" not including end sink / 10' - 9" with reagent shelf
H	36 to countertop / 54 with reagent shelf
D	54 countertop / 12 with reagent shelf

Combination options, sizes and lengths can be customized to fit your requirements.

Assemblies

Perimeter Labs




CABINET STYLE OPTIONS - All cabinets are 24"W.
The wall annex is 18"D.



201



218



242, 243



443

PERIMETER WORKSTATIONS - Provide work space for up to four people while storing your valuable laboratory equipment in these versatile perimeter workstations. Choose from four different configurations (door, door and 1 drawer, door and 4 drawers or 4 drawers). The stations are designed to utilize space along walls where they can stand alone or be ganged together. Each unit includes: 1" solid epoxy resin top with a 4" curb, solid epoxy drop-in style sink with strainer, 2 multiservice fixtures with cold water and 2 gas cocks, 4 rod sockets, rod storage, and 2 GFI receptacles. 90"W x 36"H x 84"D.

PSL-201-PL Door Unit

PSL-218-PL Door / Drawer Unit

PSL-243-PL Door / 4 Drawer Unit

PSL-443-PL 4 Drawer Unit

PERIMETER WORKSTATION WITH COUNTERTOP - Same as above except that the water and gas assemblies and the rod sockets are omitted.

PSL-201CT-PL Door Unit with Flat Top

PSL-218CT-PL Door / Drawer Unit with Flat Top

PSL-243CT-PL Door / 4 Drawer Unit with Flat Top

PSL-443CT-PL 4 Drawer Unit with Flat Top

MIX AND MATCH PERIMETER WORKSTATIONS - Provide work space for up to four people while storing your valuable laboratory equipment in these versatile perimeter workstations. Choose one style cabinet for the wall annex and another for the peninsula base. The stations are designed to utilize space along walls where they can stand alone or be ganged together. Each unit includes: 1" solid epoxy resin top with a 4" curb, solid epoxy drop-in style sink with strainer, 2 multiservice fixtures with cold water and 2 gas cocks, 4 rod sockets, rod storage, and 2 GFI receptacles. 90"W x 36"H x 84"D.

CABINET STYLE OPTIONS



201



218




242, 243



443

Assemblies



Sink Workstations



ADA SINK STATION MULTIPURPOSE - Epoxy resin countertop and sink basin. Combination hot and cold water and gas fixture. Use this item number for any ADA configuration. Add special information to drawings to meet requirements.

AS-100-PL


W	Specify
H	Specify
D	Specify



ISLAND SINK STATION - Stainless steel sinks. Two hot and two cold gooseneck faucets. Two sink cabinets with hinged doors.

AS-200-PL


W	48
H	36
D	48



All measurements in inches.

Assemblies


Sink Workstations



WALL SINK STATION - Stainless Steel basin. Hot and cold water faucet. Two sink cabinets with closure panel.

AS-300-PL


W	72
H	36" to top of sink
D	30



WALL SINK STATION - Stainless Steel basin. Three hot and cold water faucets (one ADA). Two sink cabinets with closure panel.

AS-400-PL

W	72
H	34" to top of sink
D	30




WALL SINK STATION - Stainless Steel basin. Hot and cold water faucet. Two sink cabinets with closure panel.

AS-500-PL

W	80
H	40
D	24

All measurements in inches.

Assemblies




CHALKBOARD / MARKER BOARD STORAGE - Four boards on slides. Storage units behind boards with three adjustable shelves each. Two door base cabinets, each with adjustable shelf.

AS-301-PL

W	96, 120, 144
H	84
D	25

96"W has (2) 48"W base cabinets and (2) 48" hutch cabinets
 120"W has (3) 39"W base cabinets and (3) 39" hutch cabinets
 144"W has (3) 48"W base cabinets and (3) 48" hutch cabinets



CHALKBOARD / MARKER BOARD STORAGE - Four boards on slides. Storage units behind boards with three adjustable shelves each. Two door / drawer base cabinets, each with adjustable shelf.

AS-329-PL

W	96, 120, 144
H	84
D	25

96"W has (2) 48"W base cabinets and (2) 48" hutch cabinets
 120"W has (3) 39"W base cabinets and (3) 39" hutch cabinets
 144"W has (3) 48"W base cabinets and (3) 48" hutch cabinets

All measurements in inches.